



URINE STRIPS

Urine Strips



	Kit Size	REF	Urobilinogen	Glucose	Bilirubin	Ketones	Specific Gravity	Blood	pH	Protein	Nitrite	Leukocytes	Ascorbic Acid	Microalbumin	Creatinine	CE
Urine Strip 1	100	G04001	-	✓	-	-	-	-	-	-	-	-	-	-	-	✓
Urine Strip 2	100	G04002	-	✓	-	✓	-	-	-	-	-	-	-	-	-	✓
Urine Strip 2A	100	G04002A	-	✓	-	-	-	-	-	✓	-	-	-	-	-	✓
Urine Strip 2MC	50	G11002MC	-	-	-	-	-	-	-	-	-	-	-	✓	✓	✓
Urine Strip 3A	100	G04003A	-	✓	-	-	-	-	✓	✓	-	-	-	-	-	✓
Urine Strip 4	100	G04004	-	✓	-	✓	-	-	✓	✓	-	-	-	-	-	✓
Urine Strip 4SG	100	G04004SG	-	✓	-	-	✓	-	✓	✓	-	-	-	-	-	✓
Urine Strip 5	100	G04005	-	✓	-	✓	-	✓	✓	✓	-	-	-	-	-	✓
Urine Strip 9	100	G04009	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	✓
Urine Strip 10C	100	G04010C	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	-	-	✓
Urine Strip 10ME*	100	Z10010ME	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	-	-	✓
Urine Strip 11	100	G04011	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
Urine Strip 13	100	G18013	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

* Mesh Technology: State of the art nylon mesh to prevent chemical run-over and interaction between pads

Control

	Kit Size	REF	Urobilinogen	Glucose	Bilirubin	Ketones	Specific Gravity	Blood	pH	Protein	Nitrite	Leukocytes	Ascorbic Acid	Microalbumin	Creatinine	CE
Urine Comby Control**	2x 12 mL	798001	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	✓	✓	✓

** positive and negative control for Urine Strips for following parameters

Instruments

Urine Strip Reader 40




- Throughput: 60/120 tests/h
- Closed system with barcode
- Data memory: 2000 patient results
- Built-in thermal printer
- For all urine strips except 1 parameter & 10ME
- CE marked



Urine Strip Analyser 500

- Throughput: 500 tests/h
- Large touch screen
- Closed system with barcode
- Data memory: 2000 patient results
- Built-in thermal printer
- For all urine strips except 1 parameter & 10ME
- CE marked





DIALAB®

Austria

DIALAB® GmbH

IZ NOE-Sued

Hondastrasse, Objekt M55

A-2351 Wiener Neudorf

Phone: +43-2236-66 09 10-0

Fax: +43-2236-66 09 10-30

E-mail: office@dialab.at

www.dialab.at